


## G540 OXYGEN SERIES



PATENT PENDING

FasTest's G540 gas connector is simple to use with a no twist, no turn approach to filling. Simply connect onto your CGA 540 valve and lock into place with a forward movement of the bail handle roller actuator and you are ready to fill the oxygen cylinder in a matter of seconds.

- Bail handle actuation connects to CGA 540 valve
- Available in Stainless Steel or Monel and brass. All wetted components are brass, monel, and inconel
- Cleaned for oxygen service
- Heavy duty handle eliminates common failures
- Simplified design has fewer seals to maintain
- One piece exterior body construction
- Available in rided, flexible cable, and lever style handles



Product Family	Connection Type	Connection Method	Application
G-Series	External Threads	Bail Handle	Oxygen
Connects To	Pressure Rating (psi)	Termination	Weight
CGA 540	3600	1/4" NPT (Male)	40 ounces

Part Number	Handle	Housing	Seals	Replacement Seals
G540-0411-BH	1	Stainless Steel and Brass	Viton	SG540
G540-0412-BH	2			
G540-0414-BH	4			
G540-0415-BH	5			

For EPDM seals add "E" to the part number. Example: G540-0411-BHE

G540-0411-BHM	1	Monel / Inconel and Brass	Viton	SG540
G540-0412-BHM	2			
G540-0414-BHM	4			
G540-0415-BHM	5			

For EPDM seals add "E" to the part number. Example: G540-0411-BHME

Additional configurations available, consult part number key or contact **FasTest**.