



FasTest's G580 RPV gas connector is simple to use with a no twist, no turn approach to filling. Simply connect onto your CGA 580 RPV valve and lock into place with a forward movement of the bail handle roller actuator and you are ready to fill the inert gas cylinder in a matter of seconds.

- Bail handle actuation connects to CGA 580 RPV valve
- · Heavy duty handle eliminates common failures
- Cleaned for oxygen service
- Simplified design has fewer seals to maintain
- One piece exterior body construction
- Available in ridged, flexible cable, and lever style handles







Product Family	Connection Type	Connection Method	Application
G-Series	Internal Threads	Bail Handle	Inert Gas*
Connects To	Pressure Rating (psi)	Termination	Weight
CGA 580 RPV	3600	1/4" NPT (Male)	46 ounces

<sup>\*</sup>Argon, Nitrogen, Helium or mixed gases.

Part Number	Handle	Housing	Seals	Replacement Seals
G580-0411-BRP	1	Stainless Steel/ Brass	Viton	SG580RPV
G580-0412-BRP	2			
G580-0414-BRP	4			
G580-0415-BRP	5			
G580-0416-BRP	6			

For EPDM seals add "E" to the part number. Example: G580-0411-BRPE

Additional configurations available, consult part number key or contact *FasTest*.