



FasTest's G590 gas connector is simple to use with a no twist, no turn approach to filling. Simply connect onto your CGA 590 valve and lock into place with a forward movement of the bail handle roller actuator and you are ready to fill the cylinder in a matter of seconds.

- Bail handle actuation connects to CGA 590 valve
- · Heavy duty handle eliminates common failures
- Cleaned for oxygen service
- Simplified design has fewer seals to maintain
- One piece exterior body construction
- Available in ridged, flexible cable, and lever style handles







Product Family	Connection Type	Connection Method	Application	
G-Series	Internal Threads	Bail Handle	Inert Gas	
Connects To	Pressure Rating (psi)	Termination	Weight	
CGA 590	3600	1/4" NPT (Male)	46 ounces	

Part Number	Handle	Housing	Seals	Replacement Seals
G590-0411-BH	1	Stainless Steel/ Brass	Viton	SG580

For EPDM seals add "E" to the part number. Example: G590-0411-BHE

Additional configurations available, consult part number key or contact *FasTest*.